

PERFORACIÓN COLISAS PARALELAS

Medidas	Dist Centro Ancho	Dist. Centro Largo	Coef. Perf.	Med. 2000x1000										Espesor (mm)
				0.5	0.8	1	1.5	2	2.5	3	4	5	6	
1.5 x 20	4	25	29	x	x	x								
2 x 20	5	25	31	x	x	x	x							
2 x 25	5	30	42	x	x	x	x							
2.5 x 20	6	25	33	x	x	x	x							
3 x 20	7	26	32	x	x	x	x	x						
4 x 20	8	25	40	x	x	x	x	x	x	x				
4.75 x 19.75	8.3	24.2	45	x	x	x	x	x	x	x				
5 x 20	10	25	40	x	x	x	x	x	x	x				
5 x 30	9	34	47	x	x	x	x	x	x	x				
6 x 30	11	35	47	x	x	x	x	x	x	x				
6 x 30	13	37	38	x	x	x	x	x	x	x				
7 x 20	12.5	26.5	46	x	x	x	x	x	x	x				
8 x 40	Elegir	47.5	53	x	x	x	x	x	x	x				
8 x 80	16	Elegir	Elegir	x	x	x	x	x	x	x				
10 x 40	50	50	50	x	x	x	x	x	x	x				
10 x 100	Elegir	Elegir	Elegir	x	x	x	x	x	x	x	x			
12 x 40	18	50	53	x	x	x	x	x	x	x	x			
20 x 125	Elegir	Elegir	Elegir					x	x	x	x	x	x	
20 x 200	Elegir	Elegir	Elegir					x	x	x	x	x	x	
20 x 300	Elegir	Elegir	Elegir					x	x	x	x	x	x	

